

**DURATION OF SPHERASORB VS COMPETITOR IN  
THE CATO AND CICERO 1.5 KG ABSORBER**

<b>Spherasorb Hours of use</b>	<b>Sodasorb Hours of use</b>
38.15	31.40
37.30	30.50
42.40	32.35
39.10	29.50
41.25	30.30
36.45	31.20
43.20	32.10
40.45	28.25
38.50	29.30
42.10	19.45
39.30	31.40
40.45	29.50
37.30	31.10
41.10	27.50
40.45	28.30
42.30	30.55
39.35	29.25
44.25	30.35
40.20	31.40

**Results are courtesy of the University of Dusseldorf.  
Each figure represents one absorber fill.**

**Please note, these figures are those obtained at this hospital under  
particular conditions showing that Spherasorb will last longer.  
Different conditions affect duration.**